



pneumatic instruments

for large bone, small bone and cardiothoracic surgery



**MULTIDRIVE™
MODULAR & DEDICATED HANDPIECES**



DPX-676
Sagittal saw




DPX-673
Reciprocating saw

MPX-610
Fitted with A.O. drill attachment
















Three position trigger function control. Forward, reverse and safe operating mode.



HANDPIECE SPECIFICATION

| Type code | Order code | Speed | Cannulation | Inlet connection | Operating pressure | |
|-----------|------------|----------------------|-------------|------------------|--------------------|---|
| MPX-610 | 226464 | 0 - 950 rpm | Ø 4.1 mm | MPX standard | 6 - 7 bar |  |
| | 226474 | 0 - 950 rpm | Ø 4.1 mm | Zimmer | 6 - 7 bar | |
| | 226484 | 0 - 950 rpm | Ø 4.1 mm | Aesculap | 6 - 7 bar | |
| | 226494 | 0 - 950 rpm | Ø 4.1 mm | A.O. | 6 - 7 bar | |
| DPX-676 | 226384 | 0 - 14000 cycles/min | - | - | 6 - 7 bar |  |
| DPX-673 | 226394 | 0 - 11000 cycles/min | - | - | 6 - 7 bar |  |

MPX ATTACHMENTS

| | Order code | Model | Chuck description | Final drive speed rpm | Cannulation Ø mm |
|---|------------|--------|---|-------------------------|--------------------------|
| DRILLING ATTACHMENTS | | | | | |
|  | 11710 | DX-600 | 0.4 - 6.35 mm (1/4") Jacobs | 0 - 950 | 4.1 |
| | 11720 | DX-600 | 0.4 - 4.0 mm (5/32") Jacobs | 0 - 950 | 3.1 |
|  | 11860 | DX-600 | 0.4 - 6.35 mm (1/4") Jacobs Hi-torque | 0 - 500 | 4.1 |
|  | 11740 | DX-600 | 0.4 - 3.2 mm (1/8") keyless | 0 - 950 | 2.4 |
| | 11750 | DX-600 | 0.8 - 6.35 mm (1/4") keyless | 0 - 950 | 4.1 |
|  | 11730 | DX-600 | A.O. Synthes™/ASIF - Quick release | 0 - 950 | 2.0 |
| | 11780 | DX-600 | A.O. Synthes™/ASIF Hex - Quick release | 0 - 950 | 4.1 |
| | 11770 | DX-600 | Hudson - Quick release | 0 - 950 | 4.1 |
| | 11790 | DX-600 | Trinkle - Quick release | 0 - 950 | 4.1 |
| | 12420 | DX-600 | Zimmer - Quick release | 0 - 950 | 4.1 |
| | 14480 | DX-600 | Hudson/Zimmer Combi - Quick release | 0 - 950 | 4.1 |
|  | 11760 | OX-600 | Oscillating drill - 130° oscillation arc | 0 - 950 | - |
| STANDARD REAMING ATTACHMENTS | | | | | |
|  | 11800 | RX-600 | 0 - 6.35 mm Jacobs | 0 - 260 | 4.1 |
|  | 11810 | RX-600 | Hudson | 0 - 260 | 4.1 |
| | 11820 | RX-600 | A.O. Synthes™ /ASIF/PROTEK | 0 - 260 | 4.1 |
| | 11830 | RX-600 | Zimmer | 0 - 260 | 4.1 |
| | 12900 | RX-600 | Hudson/Zimmer Combi | 0 - 260 | 4.1 |
| HIGH-TORQUE REAMING ATTACHMENTS | | | | | |
|  | 11850 | HX-600 | Hudson/Zimmer Combi | 0 - 150 | - |
| | 11840 | HX-600 | A.O. Synthes™ /ASIF/PROTEK | 0 - 150 | - |
| SAWING ATTACHMENTS FOR MPX-610 | | | | | |
|  | 12390 | CX-600 | Reciprocating saw Uses S22 type blades | 0 - 9000 cycles/min | - |
|  | 12380 | KX-600 | Large sagittal saw Uses S87 type blades | 0 - 13000 cycles/min | - |
| | 14400 | MX-600 | Mini sagittal saw - uses S83 type blades | 0 - 12000 cycles/min | - |
|  | 12300 | NX-600 | Mini sagittal saw Uses S88 type blades | 0 - 25000 cycles/min | - |
| 90° RADIOLUCENT DRILL ATTACHMENT | | | | | |
|  | 11980 | TX-600 | Accepts A.O. radiolucent drill attachment | 0 - 950 | - |
|  | | | A.O. radiolucent drill attachment Order code 17710 | | |
| PULSE LAVAGE ATTACHMENT | | | | | |
|  | 15820 | LX-620 | Accepts P20 Series pump kits | 0 - 1600 pulses/min | - |
| WIRE/PIN DRIVER ATTACHMENT | | | | | |
|  | 17250 | WX-630 | Wire driver | 0 - 950 | 0.7 - 4.0 mm capacity |

Ask for Stericut® catalogue for range of sterile packed saw blades, drill bits, rotary burs and K wires.

ERGOPLUS™ MODULAR HANDPIECES



EPV-220
Fitted with rotary bur attachment



EPV-220
Fitted with sagittal saw attachment



EPV-221
Fitted with angled wire drive attachment



The fully adjustable lever rotates, offering improved usability and ergonomics.



















EPV-220 operating lever extends for optional finger positioning.

HANDPIECE SPECIFICATION


| Control method | Type code | Order code | Speed | Operating pressure |
|----------------|-----------|------------|---------------|--------------------|
| Lever | EPV-220 | 226444 | 0 - 35000 rpm | 6 - 7 bar |
| Foot pedal | EPV-221 | 226454 | 0 - 35000 rpm | 6 - 7 bar |



EPV ATTACHMENTS

| | Order code | Model | Chuck description | Final drive speed rpm | Cannulation Ø mm |
|---|------------|--------|---|-------------------------|-----------------------|
| DRILLING ATTACHMENTS | | | | | |
|  | 14560 | DV-220 | 0 - 4 mm (5/32") Jacobs | 0 - 1300 | - |
| | 15950 | DV-225 | High speed 0.4 - 4.0 mm (5/32") Jacobs | 0 - 5000 | - |
| | 15960 | DV-225 | High speed drill A.O. | 0 - 5000 | - |
|  | 12790 | DV-220 | Offset A.O. | 0 - 1500 | 2.2 |
| | 14250 | DV-220 | Offset 0 - 4.0 mm (5/32") Jacobs | 0 - 1500 | 2.2 |
| ROTARY BUR ATTACHMENTS | | | | | |
|  | 14620 | QV-220 | Medium speed drill ISO 6360 QR | 0 - 5000 | - |
|  | 12780 | BV-220 | Accepts all 2.38 mm (3/32") burs Uses R Series burs | 0 - 35000 | - |
| | | |  Long bur guard Order code 14810  Extra long bur guard Order code 14820 | | |
|  | 17320 | BV-228 | Metal cutter - includes splash hood Accepts metal cutter R60-201 | 0 - 17500 | - |
|  | 14590 | EV-220 | E-Fitting adaptor Used with AV-220 angle rotary bur | 0 - 35000 | - |
|  | 14940 | AV-221 | 20° Contra-angle rotary bur 1:1 Used with E - Fitting adaptor 14590 | 0 - 35000 | - |
|  | 9770 | AV-220 | Standard length 20° rotary bur 1:1 | 0 - 35000 | - |
| | 9780 | AV-220 | Extended length 20° rotary bur 1:1 | 0 - 35000 | - |
| | | | 20° Contra-angle rotary bur 1:1 Used with E - Fitting adaptor 14590 | | |
| RECIPROCATING SAW ATTACHMENT | | | | | |
|  | 14550 | CV-220 | 2.5 mm stroke Uses S92 type blades | 0 - 17000 cycles/min | - |
| SAGITTAL SAW ATTACHMENTS | | | | | |
|  | 14570 | MV-220 | Micro sagittal saw - 4° oscillating arc Uses S86 type blades | 0 - 28000 cycles/min | - |
|  | 14610 | SV-220 | Sagittal saw - 4.5° oscillating arc Uses S85 type blades | 0 - 28000 cycles/min | - |
|  | 14790 | NV-220 | Sagittal saw - quick release Uses S88 type blades | 0 - 28000 cycles/min | - |
| OSCILLATING SAW ATTACHMENT | | | | | |
|  | 14580 | OV-220 | Wrenchless blade chucking Uses S86 type blades | 0 - 22000 cycles/min | - |
| WIRE DRIVE ATTACHMENT | | | | | |
|  | 17010 | WV-220 | Angled wire driver | 0 - 1200 | 0.6 - 1.6 capacity |

Ask for Stericut® catalogue for range of sterile packed saw blades, drill bits, rotary burs and K wires.



MPZ-450
Fitted with sagittal saw attachment


MPZ-450
Fitted with quick release drill attachment



MPZ-451 for foot pedal operation



*MPZ handpiece wash cap
Order code 9950*





















See page 9 for foot pedal specifications

HANDPIECE SPECIFICATION

| Control method | Type code | Order code | Speed | Operating pressure |
|----------------|-----------|------------|--------------|--------------------|
| Trigger | MPZ-450 | 226414 | 0 - 1100 rpm | 6 - 7 bar |
| Foot pedal | MPZ-451 | 226424 | 0 - 1100 rpm | 6 - 7 bar |



MPZ ATTACHMENTS

| | Order code | Model | Chuck description | Final drive speed rpm | Cannulation Ø mm |
|---|------------|---|--|--|------------------|
| DRILLING ATTACHMENTS | | | | | |
|  | 14110 | DZ-450 | 0 - 4.0 mm (5/32") Jacobs | 0 - 1100 | 3.1 |
| | 14240 | DZ-450 | 0 - 6.35 mm (1/4") Jacobs | 0 - 1100 | 4.1 |
|  | 14130 | DZ-450 | 0 - 3.0 mm (1/8") Keyless | 0 - 1100 | 3.2 |
| | 14120 | DZ-450 | 0 - 6.35 mm (1/4") Keyless | 0 - 1100 | 4.1 |
| QUICK RELEASE DRILL ATTACHMENTS | | | | | |
|  | 14140 | DZ-450 | A.O. Synthes™/ASIF | 0 - 1100 | 2.0 |
| | 15040 | DZ-450 | Hudson/Zimmer Combi | 0 - 1100 | 4.1 |
| | 15050 | DZ-450 | Trinkle | 0 - 1100 | 4.1 |
|  | 15800 | QZ-450 | Medium speed drill ISO 6360 QR | 0 - 7000 | - |
| HIGH SPEED BUR | | | | | |
|  | 16900 | BZ-450 | Accepts all 2.38 mm (3/32") burs Uses R Series burs | 0 - 22000 | - |
| | | |  Long bur guard Order code 14810 |  Extra long bur guard Order code 14820 | |
| REAMING ATTACHMENTS | | | | | |
|  | 14160 | RZ-450 | 0 - 6.35 mm (1/4") Jacobs | 0 - 300 | 4.1 |
|  | 14150 | RZ-450 | A.O. Synthes™/ASIF/PROTEK | 0 - 300 | 4.1 |
| | 14170 | RZ-450 | Hudson/Zimmer | 0 - 300 | 4.1 |
| SAWING ATTACHMENTS | | | | | |
|  | 14210 | SZ-450 | Sagittal saw Uses S85 type blades | 0 - 18000 cycles/min | - |
|  | 14230 | NZ-450 | Sagittal saw - quick release Uses S88 type blades | 0 - 21000 cycles/min | - |
|  | 14220 | CZ-450 | Reciprocating saw Uses S92 type blades | 0 - 8000 cycles/min | - |
|  | 15110 | CZ-451 | Reciprocating saw Uses S22 type blades | 0 - 6000 cycles/min | - |
| RADIO LUCENT ADAPTOR | | | | | |
|  | 15210 | TZ-450 | Accepts A.O. radiolucent drill attachment | 0 - 1100 | - |
| PULSE LAVAGE ADAPTOR | | | | | |
|  | 17230 | LZ-451 | Pulse lavage attachment Uses P20-7XX Series pump kits | 0 - 1100 pulses/min | - |
| WIRE DRIVE ATTACHMENT | | | | | |
|  | 17190 | WZ-451 | Wire/pin driver | 0 - 1100 | 0.7 - 4.0 |
| | 14260 | WZ-450 | Wire driver - gold nose | 0 - 1100 | 0.6 - 1.8 |
| | 14270 | PZ-450 | Pin driver - silver nose | 0 - 1100 | 1.8 - 3.2 |
| | |  | | | |
| | |  | | | |

Ask for Stericut® catalogue for range of sterile packed saw blades, drill bits, rotary burs and K wires.

MICRODRIVE DEDICATED HANDPIECES



D30-100
Drill/bur



D32-100
Sagittal saw



D33-100
Oscillating saw



D34-100
Reciprocating saw



D36-100
Wire driver



Extra long bur guard
Order code 2240

Long bur guard
Order code 2230



Handpiece wash cap
Order code 8110

HANDPIECE SPECIFICATION

| Function | Type code | Order code | Speed | Wire capacity/ cannulation | Air supply operating pressure | Accessories |
|-------------------|-----------|------------|----------------------|-------------------------------|----------------------------------|-----------------------------|
| Drill | D30-100 | 219664 | 0 - 30000 rpm | - | 6 - 7 bar | Stericut R Series burs |
| Sagittal saw | D32-100 | 219774 | 0 - 25000 cycles/min | - | 6 - 7 bar | Stericut S92 saw blades |
| Oscillating saw | D33-100 | 219884 | 0 - 25000 cycles/min | - | 6 - 7 bar | Stericut S92 saw blades |
| Reciprocating saw | D34-100 | 219994 | 0 - 18000 cycles/min | - | 6 - 7 bar | Stericut S92 saw blades |
| Wire driver | D36-100 | 226004 | 0 - 2000 rpm | 0.7 - 1.6 mm | 6 - 7 bar | Stericut K wires and guides |

PNEUMATIC INSTRUMENT ACCESSORIES

IPV-220 IRRIGATION PUMP



Order code 11560


Used with EPV-220 or EPV-221 via air foot pedal.

- Fully disposable irrigation kits.
- Simple 'clip-on set-up.
- 4 function settings.
- Variable flow rate.
- Removable pump head for cleaning.

Input - 240Vac, 120Vac and 110Vac

Mains power cord sets - 3 metres long. Please specify

| | |
|-----------------------|------------------|
| UK | Order code 8290 |
| Euro | Order code 8300 |
| US/Canada | Order code 8310 |
| Australia/New Zealand | Order code 8320 |
| Swiss | Order code 14700 |



Irrigation hose kit for use with rotary bur applications
Order code T10-103

PUMP SENSOR AND AIR HOSE FOR IPV-220



MPZ Wash Assist air hose sub assembly

Hose length - Order code

| 3 metre | 4 metre | 5 metre |
|---------|---------|---------|
| 632333 | 632343 | 632353 |




IPS-220 Pump sensor

Order code 11570

Detects small amounts of air flow in the exhaust to start the pump. Fits simply between the air hose connecting to the EPV-221 and the air supply or foot pedal.

AIR FOOT PEDAL



| Model | Order code | Foot pedal features | |
|--------|------------|--|---|
| | | Inlet connection | Outlet connection |
| FV-220 | 14900 | Integral 3 metre supply hose with MA7 probe | Accepts any MPZ "Wash Assist" air hose sub assembly - see above |
| FV-221 | 14910 | Accepts Draeger Wash Assist hoses with air being exhausted at wall fitting | Accepts any MPZ "Wash Assist" air hose sub assembly - see above |
| FV-222 | 15000 | Integral 3 metre supply hose with MA7 probe | MA7 outlet |

AIR STERNUM SAWS

DPC-470 offers the same innovative functionality as the battery sternum saw. A single piece, twist and click, micro adjustable blade guard provides quick and efficient use with optimum control of the cutting direction. For revision sternotomies, the DPX-170 oscillating saw handpiece is also offered.



DPC-470
for primary sternotomies



*Simple twist and click
sternum blade guard for
quick and efficient use.*

HANDPIECE SPECIFICATION

| Model | Order code | Speed | Stroke | Operating pressure | Sternum guard Order code |
|---------|------------|---------------|--------|--------------------|--------------------------|
| DPC-470 | 226614 | 0 - 16000 cpm | 3.7 mm | 6 - 7 bar | 16490 |



| Model | Order code | Speed | Oscillating arc | Operating pressure |
|---------|------------|---------------|-----------------|--------------------|
| DPX-170 | 223224 | 0 - 16000 cpm | 6° | 6 - 7 bar |



Air inlet washing cap.
Prevents water entering handpiece
during washing prior to autoclaving.



Order code 9950 for DPC-470
Order code 608073 for DPX-170


ECOPULSE™ PULSE LAVAGE SYSTEM



PLX-620
Fitted with P20 pump kit

ECO Pulse pump kits simply “click” into any PLX rotary handpiece

Up to 60% less clinical waste, helping to conserve valuable hospital storage by up to 50%.



| NOZZLES | | | |
|--|--------|------------|--|
| Type | Length | Order code | |
| Soft tissue less suction | 10 cm | P20-120 | |
| Soft tissue with suction | 10 cm | P20-121 | |
| Bone cleaning with suction | 10 cm | P20-125 | |
| Femoral canal less suction | 25 cm | P20-122 | |
| Femoral canal with suction | 25 cm | P20-123 | |
| Femoral canal (high capacity) with suction | 25 cm | P20-126 | |
| Femoral brush with suction | 25 cm | P20-124 | |

All pump kits are supplied with nozzles and 3 metres of tubing.
Suction tubing is supplied with all suction types.
All sterile packed pump kits and individual nozzles are supplied 5 per box.

| PUMP KITS | | Type | Order code | Supplied nozzles |
|-----------|--|--|------------|------------------|
| | | General purpose with suction | P20-100 | P20-121 |
| | | Hip arthroplasty with suction | P20-101 | P20-121/P20-123 |
| | | Hip arthroplasty brush with suction | P20-102 | P20-124 |
| | | General purpose without suction | P20-103 | P20-120 |
| | | Hip arthroplasty without suction | P20-104 | P20-120/P20-122 |
| | | Bone cleaning (high pressure) with suction | P20-105 | P20-125 |
| | | Bone cleaning (high pressure) without suction | P20-106 | P20-125 |
| | | Hip arthroplasty (high pressure) with suction | P20-107 | P20-125/P20-123 |
| | | Hip arthroplasty (high pressure) without suction | P20-708 | P20-125/P20-122 |

HANDPIECE SPECIFICATION

| Model | Order code | Air inlet | Variable pulse rate Pulses/min | Variable flow rate Litres/min | Operating pressure bar |
|---------|------------|--------------|-----------------------------------|----------------------------------|---------------------------|
| PLX-620 | 226574 | PLX standard | 0 - 1450 | 0 - 1.20 | 6 - 7 |
| PLX-621 | 226584 | MPX | 0 - 1450 | 0 - 1.20 | 6 - 7 |
| PLX-623 | 226604 | A.O. | 0 - 1450 | 0 - 1.20 | 6 - 7 |



STANDARD CONTAINER AND UNIVERSAL TRAY



Universal sterilisation case and tray

Order code: 9350


Universal tray Order code: 631343

Lid Order code: 630843


Base Order code: 630833

Blue mat Order code: 616153

Will accept all
standard trays



WIRE CONTAINER




Wire sterilisation container and lid.

Order code 655073

Accepts both sterile and aseptic insert trays

STERILE FILTERED CONTAINER



Order code 17410

Accepts standard or wire container plus tray



Filter paper - single use

pack of 100

Order code 657243




Security tag,

pack of 100

Order code 657253

Sterile indicator label, pack of 100 Order code 657263

HALF SIZE WIRE CONTAINER




Wire Sterilisation Container and lid.

Order code 655593

STERILISATION EQUIPMENT

MULTIDRIVE™ PNEUMATIC SYSTEM



| Set up | Standard Container and insert tray Order code | Insert tray Order code |
|----------------------|---|------------------------|
| 3 MPX/DPX handpieces | 14630 | 14660 |
| 2 MPX/DPX handpieces | 14640 | 14650 |



Washing cap for MPX, SPX and DPX handpieces
Order code 608073

MULTIDRIVE™ PNEUMATIC SYSTEM



| Set up | Standard Container and insert tray Order code | Insert tray Order code |
|---|---|------------------------|
| 1 MPZ and 1 EPV handpiece, 5 EPV attachments, 8 MPZ attachments and an air hose | 14390 | 14410 |



Washing cap for EPV, MPZ, HST and DPZ handpieces
Order code 8110

STERNUDRIVE™ PNEUMATIC SYSTEM



| Set up | Standard Container and insert tray Order code | Insert tray Order code |
|---------|---|------------------------|
| DPC-470 | 17160 | 17170 |
| DPX-170 | 9350 * | 631343 |

Insert tray silicone mat - Order code 616153
* Refer to page 10



Washing cap for DPX-170 handpiece
Order code 608073



Washing cap for DPC-470 handpiece
Order code 9950

MICRODRIVE™ PNEUMATIC SYSTEM



D30 Series handpieces
Order code 7540











Washing cap for D30 Series and PLX handpieces
Order code 8110

Insert trays and Containers shown complete with handpieces and other items are for illustration purposes only.

WASH ASSIST AIR HOSES

The Wash Assist co-axial hose system is recommended where machine washing of air hoses is demanded. By releasing the probe end from hose, which is retained on the hose via the rubber strap, the hose ends can be joined and sealed prior to washing.








| Probe end | For model type | Complete hose length - Order code | | |
|---|--|-----------------------------------|---------|---------|
| | | 3 metre | 4 metre | 5 metre |
|  MA7-Schrader | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13000 | 13100 | 13200 |
| | | 13300 | 13400 | 13500 |
| | | 13600 | 13700 | 13800 |
|  A.O. | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13010 | 13110 | 13210 |
| | | 13310 | 13410 | 13510 |
| | | 13610 | 13710 | 13810 |
|  AGA Inst. | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13020 | 13120 | 13220 |
| | | 13320 | 13420 | 13520 |
| | | 13620 | 13720 | 13820 |
|  Draeger | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13030 | 13130 | 13230 |
| | | 13330 | 13430 | 13530 |
| | | 13630 | 13730 | 13830 |
|  Heyer | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13040 | 13140 | 13240 |
| | | 13340 | 13440 | 13540 |
| | | 13640 | 13740 | 13840 |
|  Air Liquide | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13050 | 13150 | 13250 |
| | | 13350 | 13450 | 13550 |
| | | 13650 | 13750 | 13850 |
|  DIN | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13060 | 13160 | 13260 |
| | | 13360 | 13460 | 13560 |
| | | 13660 | 13760 | 13860 |
|  AGA Air | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13080 | 13180 | 13280 |
| | | 13380 | 13480 | 13580 |
| | | 13680 | 13780 | 13880 |
| Flowtech | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13070 | 13170 | 13270 |
| | | 13370 | 13470 | 13570 |
| | | 13670 | 13770 | 13870 |
| NP05 | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 13090 | 13190 | 13290 |
| | | 13390 | 13490 | 13590 |
| | | 13690 | 13790 | 13890 |

To special order

STANDARD AIR HOSES



Quick release, bayonet action hose fitting with combined coaxial air supply and exhaust hose. This system is non-preferred for machine washing.

| Probe end | For model type | Complete hose length - Order code | | |
|---|--|-----------------------------------|---------|---------|
| | | 3 metre | 4 metre | 5 metre |
|  | MA7-Schrader MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 5300 | 5310 | 5320 |
| | | 10900 | 10910 | 10920 |
| | | 7500 | 9320 | 7640 |
|  | A.O. MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 5330 | 5340 | 5350 |
| | | 10930 | 10940 | 10950 |
| | | 8690 | 8700 | 8710 |
|  | AGA Inst. MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 5360 | 5370 | 5380 |
| | | 10960 | 10970 | 10980 |
| | | 8720 | 8730 | 8740 |
|  | Draeger MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 5390 | 5400 | 5410 |
| | | 10990 | 11000 | 11010 |
| | | 8580 | 8590 | 8600 |
|  | Heyer MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series | 7820 | 7830 | 7840 |
| | | 11020 | 11030 | 11040 |
| | | 8750 | 8760 | 8770 |
| Air Liquide | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 | 10390 | 10400 | 10410 |
| | | 11050 | 11060 | 11070 |
| Flowtech | MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 | 8610 | 8620 | 8630 |
| | | 11080 | 11090 | 11100 |

To special order

HOSE WASHING CAP



| Probe end | For model type | | |
|--------------|------------------------------|---|-------------------------------|
| | MPX, DPX, SPX, Order code | MPZ, EPV, HST, DPZ-121, DPC-470 Order code | PLX, D30 Series Order code |
| MA7 Schrader | 8950 | 9650 | 11430 |
| A.O. | 10880 | 9670 | 11450 |
| AGA | 10870 | 9660 | 11440 |
| Draeger | 10890 | 9680 | 11460 |
| Heyer | 10870 | 9660 | 11440 |

UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)
De Soutter Medical Limited
Halton Brook Business Park, Weston Road
Aston Clinton, Aylesbury, Bucks HP22 5WF
Tel: +44 (0) 1296 634 000
Fax: +44 (0) 1296 634 033
Email: info@de-soutter.com

AUSTRALIA
De Soutter Medical Australia Pty Ltd
2/12-14 Apollo Drive, Hallam
Victoria 3803
Tel: +61 (0) 3 9702 4441
Fax: +61 (0) 3 9702 4484
Email: australia@de-soutter.com

AUSTRIA
De Soutter Medical Austria
Zweigniederlassung Österreich
Dietrichsteingasse 10, A-3400 Klosterneuburg
Tel: +43 (0) 676 96 71 770
Fax: +43 (0) 2243 21 656
Email: austria@de-soutter.com

BELGIUM
De Soutter Medical Belgium
In De Bruel 30
3620 Lanaken
Tel: +32 (0) 89/47 15 37
Fax: +32 (0) 89/70 12 19
Email: belgium@de-soutter.com

FRANCE
De Soutter Medical France
949 Avenue du Parc des Expositions
33260 La Teste de Buch
Tel: +33 (0) 5 56 54 89 36
Fax: +33 (0) 9 70 61 37 60
Email: france@de-soutter.com

GERMANY
De Soutter Medical Germany
Niederlassung Deutschland, Kladenfloss
66625 Nohfelden
Tel: +49 (0) 68 52-99 12 46
Fax: +49 (0) 68 52-99 12 47
Email: deutschland@de-soutter.com

ITALY
De Soutter Medical Italia
Località Fornace SNC
27022 Casorate Primo - PV
Tel: +39 (0) 2 9009 4098
Fax: +39 (0) 2 9009 2673
Email: italy@de-soutter.com

NETHERLANDS
De Soutter Medical Netherlands
Gelderlandhaven 2X
3433 PG Nieuwegein
Nederland
Email: netherlands@de-soutter.com

UNITED STATES OF AMERICA
De Soutter Medical USA, Inc.
224 Rolling Hill Road, Suite 12A
Mooresville, NC 28117
Tel: +1 (704) 655 9040
fax: +1 (704) 987 2035
Email: usa@de-soutter.com

deSoutter
MEDICAL

For further information on De Soutter Medical powered instruments and cast saw systems please contact the address in your area.